

Product Description

Pioneering GTPase and Oncogene Product Development since 2010

EGFR(G863D)

Catalog Number: 26121

Gene Symbol: EGFR/HER1/ERBB1

Description: Anti-EGFR(G863D) Mouse Monoclonal Antibody

Background: EGFR is a receptor tyrosine kinase binding ligands of the EGF family and activating several signaling cascades to convert extracellular cues into appropriate cellular responses. It activates at least 4 major downstream signaling cascades including the RAS-RAF-MEK-ERK, PI3 kinase-AKT, PLCgamma-PKC and STATs modules. Defects in EGFR are associated with lung cancer (LNCR).

Immunogen: A synthetic peptide from the internal region of EGFR which includes the

mutation of G863D, human origin. **Tested Applications:** ELISA, WB, IHC

Recommended Dilutions:

ELISA: 1:1000-1:5000 WB: 1:500-1:1000 IHC: 1:50-1:100

Concentration: 1 mg/ml

Host: Mouse

Clonality: Monoclonal

Isotype: IgG

Purity: Purified from ascites

Format: Liquid Preservative: No

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH7.4, 150


mM NaCl, 50% glycerol

Species Reactivity: Recognizes EGFR(G863D), but

not wild type EGFR from vertebrates. **Storage Conditions:** Store at -20°C. Avoid

repeated freezing and thawing

Western blot:



1 2 WB: Anti-His

Western blot analysis of recombinant EGFR (G863D) and wildtype proteins.


Purified His-tagged EGFR(G863D) protein (amino acid702-987) (Cat. # 10222) and corresponding wild-type protein (Cat. #10210) were blotted withanti-EGFR (G863D) mouse monoclonal antibody(Cat. # 26119). Proteins loading were blotted with anti-His-tag monoclonal antibody (Cat. # 28004).

info@neweastbio.com Web: www.neweastbio.com



Pioneering GTPase and Oncogene Product Development since 2010

Immunofluorescence:



Immunofluorescence of cells expressing EGFR proteins with anti-EGFR (G863D) antibody.

HEK293T cells were transfected with pCDNA3-GFP-EGFR (WT) plasmid (left column) or pCDNA3-GFP-EGFR (G863D) plasmid (right column), then fixed and stained with anti-EGFR (G863D) monoclonal antibody (Cat. # 26119).

www.neweastbio.com