

Product Description


Pioneering GTPase and Oncogene Product Development since 2010

ALK(R1275Q)

Catalog Number: 26151 Gene Symbol: ALK, CD246, NBLST3 Description: Anti-ALK(R1275Q) Mouse Monoclonal Antibody

Background: Anaplastic lymphoma kinase (ALK) also known as CD246 (cluster of differentiation 246) is a tyrosine kinase receptorthat plays an important role in the development of the brain and nervous systems. Deficiencies in ALK, which most commonly as a fusion withother genes, are indicated in cancers. Immunogen: A synthetic peptide from theinternal region of ALK, which includes the mutation of R1275Q, human origin. Tested Applications: ELISA, WB, IHC **Recommended Dilutions:** ELISA: 1:1000-1:5000 WB: 1:500-1:2000 IHC: 1:50-1:100 Concentration: 1 mg/ml Host: Mouse Clonality: Monoclonal Isotype: IgG Purity: Purified from ascites Format: Liquid Preservative: No Constituents: PBS (without Mg²⁺ and Ca²⁺), pH7.4, 150 mM NaCl, 50% glycerol Species Reactivity: Recognizes R1275Q mutant, but not wild type ALK of vertebrates. Storage Conditions: Store at -20°C. Avoid repeated freezing and thawing

Western blot:



Western blot analysis of recombinant ALK (R1275Q) and wild type proteins.

Purified His-tagged ALK (R1275Q) protein (amino aci1060-1620, lane 2) and corresponding wild type protein (lane 1) were blotted with anti-ALK (R1275Q) mouse monoclonal antibody (Cat. # 26151).

FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC APPLICATIONS

NewEast Biosciences 1150 First Avenue, Suite 501 King of Prussia, PA 19406 Support: Web:

610-945-2007 info@neweastbio.com www.neweastbio.com