

Ras(G12S)

Catalog Number: 26186

Gene Symbol: H-Ras; K-Ras; N-Ras

Description: Anti-RAS (G12S) Mouse Monoclonal Antibody

Background: Ras protein is a membrane-associated GTPase signaling protein that regulates proliferation, differentiation, and cell survival. Ras G12S mutation result in decreased GTPase activity and constitutive signaling and is correlated with the development of pancreatic intraepithelial neoplasia and lung lesions.

Immunogen: A synthetic peptide from the internal region of RAS which includes the mutation of G12S, human origin.

Tested Applications: ELISA, WB, IHC

Recommended Dilutions:

ELISA: 1:1000-1:5000

WB: 1:500-1:1000

IHC: 1:50-1:100

Concentration: 1 mg/ml

Host: Mouse

Clonality: Monoclonal

Purity: Purified from ascites

Format: Liquid


Preservative: No

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH7.4, 150 mM NaCl, 50% glycerol

Species Reactivity: Recognizes G12S mutant, but not wild-type RAS of vertebrates.

Storage Conditions: Store at -20°C. Avoid freeze / thaw cycles.

Western blot:




Western blot analysis of recombinant RAS (G12S) and wildtype proteins.

Purified His-tagged RAS(G12S) protein (lane 2) and corresponding wildtype protein (lane 1) were blotted with anti-RAS(G12S) monoclonal antibody (Cat. # 26186).

FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC APPLICATIONS

Immunofluorescence:



Immunofluorescence of cells expressing RAS proteins with anti-RAS(G12S) antibody.

HEK293T cells were transfected with pCDNA3-GFP-RAS (WT) plasmid (left column) or pCDNA3-GFP-RAS(G12S) plasmid (right column), then fixed and stained with anti-RAS(G12S) monoclonal antibody (Cat. # 26186).

FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC APPLICATIONS