

EGFR(L858Q)

Catalog Number: 26221

Gene Symbol: EGFR/HER1/ERBB1

Description: Anti-EGFR(L858Q) Mouse Monoclonal Antibody

Background: EGFR is a receptor tyrosine kinase binding ligands of the EGF family and activating several signaling cascades to convert extracellular cues into appropriate cellular responses. It activates at least 4 major downstream signaling cascades including the RAS-RAF-MEK-ERK, PI3 kinase-AKT, PLCgamma-PKC and STATs modules. Defects in EGFR are associated with lung cancer (LNCR).

Immunogen: A synthetic peptide from the internal region of EGFR, which includes the mutation of L858Q, human origin.

Tested Applications: ELISA, WB, IF, IHC

Recommended Dilutions:

ELISA: 1:1000-1:5000

WB: 1:500-1:1000

IHC: 1:50-1:100

Concentration: 1 mg/ml

Host: Mouse

Clonality: Monoclonal

Purity: Purified from ascites

Format: Liquid


Preservative: No

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.4, 150 mM NaCl, 50% glycerol

Species Reactivity: Recognizes EGFR (L858Q), but not wild type EGFR protein from vertebrates.

Storage Conditions: Store at -20°C. Avoid repeated freezing and thawing

Western blot:



WB: Anti-EGFR(L858Q) mAb

Western blot analysis of recombinant EGFR(L858Q) and wild type proteins.

Purified His-tagged EGFR (L858Q) protein (aa702-987, lane 1) and corresponding wild type protein (aa702-987, lane 2) were blotted with anti-EGFR(L858Q) mouse antibody (Cat. # 26221).

FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC APPLICATIONS