

Product Description

A helping hand for your research

EML4-ALK

Catalog Number: 26226 Gene Symbol: EML4-ALK

Description: Anti-EML4-ALK Mouse Monoclonal Antibody

Background: EML4-ALK is a fusion-type protein kinase Western blot: that can be found in many tumors, such as anaplastic large cell lymphoma, inflammatory myofibroblastic tumor, neuroblastoma, and NSCLC, and so on. It is generated as the result of a small inversion within the short arm of human chromosome 2.

EML4-ALK protein often promotes and maintains the malignant behavior of the cancer cells by activating the MAPK, PI3K/AKT and JAK/STAT3 pathways.

Immunogen: A synthetic peptide from the internal

region of EML4-ALK of human origin Tested Applications: ELISA, WB, IF, IHC

Recommended Dilutions:

1:1000-1:5000 ELISA: WB: 1:500-1:2000 IF: 1:50-1:100 IHC: 1:50-1:100

Concentration: 1 mg/ml

Host: Mouse

Clonality: Monoclonal

Isotype: IgG

Purity: Purified from ascites

Format: Liquid Preservative: No

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH7.4, 150


mM NaCl, 50% glycerol

Species Reactivity: Recognizes T41A mutant, but not

wild-type EML4-ALK of vertebrates.

Storage Conditions: Store at -20°C. Avoid freeze / thaw

cycles.



WB: Anti- EML4-ALK mAb Western blot analysis of recombinant EML4-ALK and wild type proteins.

PurifiedHis-tagged EML4-ALK protein (lane 3), corresponding wild-type EML4 protein (lane1) and wild-type ALK protein (lane2) were blotted with anti-EML4-ALK monoclonal antibody (Cat. # 26242).

NewEast Biosciences 1150 First Avenue, Suite 501 King of Prussia, PA 19406


610-945-2007 Support:

info@neweastbio.com Web: www.neweastbio.com



A helping hand for your research

Immunofluorescence:



Immunofluorescence of cells expressing EML4-ALK proteins with anti-EML4-ALK antibody.

HEK293T cells were transfected with pCDNA3-GFP-EML4 (WT) plasmid (left column) or pCDNA3-GFP-ALK (WT) plasmid (middle column) or pCDNA3-GFP-EML4-ALK plasmid (right column), then fixed and stained with anti-EML4-ALK monoclonal antibody (Cat. #26242).

Web: <u>info@neweastbio.com</u> www.neweastbio.com